# GLDH VET test kit

## For veterinary use only!

Veterinary test kit for quantitative in vitro determination of GLDH in serum or lithium heparin plasma with a solo analyser.

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**Test kit preparation: allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Return the test kit package to the refrigerator.**

## Summary

Glutamatdehydrogenase (GLDH) is a liver specific enzyme which is located in mitochondria of centrolobular hepatocytes. Due to this localization, GLDH dominantly increases in hypoxic states of the liver as well as in pronounced necrosis of hepatocytes.

Increased release of GLDH from hepatocytes can further be observed in i.e. inflammatory, neoplastic, or infectious liver pathology as well as in all diseases leading to liver cell degeneration.

## Method

Glutamatdehydrogenase catalyses the reduction of oxolutarat and at the same time the oxidation of NADH to NAD.

The reaction speed of the absorption increase is directly proportional to the GLDH activity which is measured on the SOLO system at 340 nm.

## Measurement Range

6 - 120 U/L (Lot dependent)

## Sample Material

Use 100 µl lithium heparin plasma or serum samples. (Avoid lipaemic samples)

In case samples need to be kept, store in a dark place only!

Sample material stability:

Stored at 20 °C:	4 hours
Stored at 2 - 8 °C:	1 day
Stored at - 20 °C:	4 weeks

## Test Kit

ERS cuvette filled with: Trietanolamin, alpha ketoglutarat, ammoniumacetat EDTA, ADP, LDH.  
ERS cap filled with: NADH

## Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 - 8 °C.  
Opening the vacuum package may limit the reagent stability to three months (stored at 2 - 8 °C) from the date of opening. DO NOT FREEZE!

## Warnings and Precautions

DO NOT INGEST! Avoid contact with skin and eyes. Observe all necessary precautions for the use of laboratory reagents.

## Waste Management

Always refer to local legal requirements.

## Reference Range

Dog:	< 13 U/L
Cat:	< 10 U/L
Rabbit:	< 19 U/L
Guinea pig:	< 17 U/L
Horse:	< 13 U/L
Bovine:	< 15 U/L
Sheep:	< 11 U/L
Pig:	< 7 U/L
Goat:	< 20 U/L

It is recommended that each laboratory establishes its own reference ranges.

## Quality Control

For internal quality control the Eurolyser GLDH VET control kit is recommended. Order number: VT 2000

## Precision

Reproducibility day-to-day:  
Control; N = 20; mean = 24.4 U/L; CV = 8.02%;

## Correlation

y (solo GLDH) = 1.0414 x (COBAS GLDH) + 4.2329;  
R<sup>2</sup> = 0.9494.

## Interferences

The test system has been analysed for various interferences. Criterion was the recovery within 10% of initial values.

Ascorbic Acid	30 mg/dl
Bilirubin	60 mg/dl
Haemoglobin	500 mg/dl
Triglyceride	400 mg/dl

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---

### **ATTENTION!**

Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C)!

### **1. Preparing the test system**

- 1.1 Insert RFID card
- 1.2 Place ERS cuvette into test kit rack
- 1.3 Place ERS cap into test kit rack
- 1.4 Press "Measurement" button and enter all required information using the touchscreen

---

### **2. Sample preparation**

- 2.1 Aspirate 100 µl sample material from centrifuged sample tube

---

### **3. Test processing**

- 3.1 Dispense 100 µl sample **INTO THE LIQUID** in the ERS cuvette
- 3.2 Firmly apply ERS cap on ERS cuvette
- 3.3 Place ERS cartridge into the solo analyser
- 3.4 Start automatic sample processing by pressing the  start button on the solo analyser