Lipase (pancreatic) VET test kit

For veterinary use only!

Veterinary test kit for quantitative in vitro determination of Lipase (pancreatic) from serum or lithium heparin plasma with a solo analyser.

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Kit size

16 tests/kit
6 tests/kit
1 x 5 ml

Test kit preparation: allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Return the test kit package to the refrigerator.

Summary

This assay uses the lipase enzymatic method DGGR which is well known as a high sensitive and specific method for the detection of pancreatitis in animals. Using this assay, Lipase is a highly effective biomarker and has shown a > 95% concordance by comparison with a tissue- and species-specific qualitative immunoassay (O'Brien et al., proceedings 15th annual congress of the ESVCP, 2013, p. 125).

Method

In the presence of colipase and bile acids lipase splits the synthetic substrate (1,2-o-dilauryl-rac-glycero-3-glutaric acid-(6-methylresorufin) ester) to glycerol and methylresorufin-ester, which is spontaneously degraded to glutaric acid and methylresorufin. The combination of colipase and bile acid makes the reaction specific for pancreatic lipase without interference of esterases and lipolytic enzymes.

The measured absorbance is proportional to the lipase activity in the sample. Photometric measurement of the absorbance rate (kinetic) at 546 nm wavelength.

Measurement Range

When using a 20 µl pipette: 25 - 300 U/L (Lot dependent)
When using a 5 µl pipette: 250 - 1200 U/L (Lot dependent)

Sensitivity: 25 U/L

Sample material

Serum or lithium-heparin plasma. Sample volume according to the menu settings on the laboratory photometer and according to the use of a 20 µl or 5 µl pipette.

Store samples protected from light!

Sample material stability:

Stored at 2 - 8 °C: 2 days

Stored at - 20 °C: 1 year

Test Kit

ERS cuvette filled with Good's buffer.

ERS cap filled with tartrate buffer, taurodehydroxycholate and color substrate.

Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 - 8 °C.

Opening the vacuum package may limit the reagent stability to three months (stored at 2 - 8 °C) from the date of opening.

DO NOT FREEZE!

Warnings and Precautions

DO NOT INGEST! Avoid contact with skin and eyes. Observe all necessary precautions for the use of laboratory reagents.

Waste Management

Always refer to local legal requirements.

Reference Range

Cat: 0 - 35 U/L

Dog: 0 - 125 U/L

It is recommended that each laboratory establishes its own reference ranges.

Quality Control

For internal quality control the Eurolyser Lipase VET control kit is recommended. Order number: VT 1800

Precision

Reproducibility within-run:

Dog serum; N = 20; mean = 137 U/L; CV = 3%;

Correlation

Dog sample correlation:

y (solo Lipase) = 1.2286 x (COBAS Lipase) - 9.0493; R² = 0.9668;

Cat sample correlation:

y (solo Lipase) = 1.1526 x (COBAS Lipase) + 4.8579; R² = 0.9713;

Interferences

The test system has been analysed for various interferences. Criterion was the recovery within 10% of initial values.

Bilirubin 60 mg/dl

Haemoglobin 500 mg/dl

Triglycerides 1000 mg/dl

References

1. Kook et al., Agreement of Serum Spec cPL with the 1,2-o-Dilauryl-Rac-Glycero Glutaric Acid-(6' -methylresorufin) Ester (DGGR) Lipase Assay and with Pancreatic Ultrasonography in Dogs with Suspected Pancreatitis, J Vet Intern Med. 2014;28:863-870;
2. Oppliger et al., Agreement of the Serum Spec cPL TM and 1,2-O-Dilauryl-Rac-Glycero-3-Glutaric Acid-(6' -Methylresorufin) Ester Lipase Assay for the Determination of Serum Lipase in Cats with Suspicion of Pancreatitis, J Vet Intern Med. 2013;27:1077-1082;
3. Gruber, Evaluation Report - Eurolyser LIPASE test kit (VT0190) on solo Analyser, Eurolyser Diagnostica GmbH, 03-06/2014;

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Bilirubin 60 mg/dl
Hämoglobin 500 mg/dl
Triglyzeride 1000 mg/dl


Referenzen

- Kook et al., Agreement of Serum Spec cPL with the 1,2-o-Dilauryl-Rac-Glycero Glutaric Acid-(6' -methylresorufin) Ester (DGGR) Lipase Assay and with Pancreatic Ultrasonography in Dogs with Suspected Pancreatitis, J Vet Intern Med. 2014;28:863-870;
- Oppliger et al., Agreement of the Serum SpecPL TM and 1,2-O-Dilauryl-RacGlycero-3-Glutaric Acid-(6' -Methylresorufin) Ester Lipase Assay for the Determination of Serum Lipase in Cats with Suspicion of Pancreatitis, J Vet Intern Med. 2013;27:1077-1082;
- Gruber, Evaluation Report - Eurolyser LIPASE test kit (VT0190) on solo Analyser, Eurolyser Diagnostica GmbH, 03-06/2014;

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Performing a Lipase (pancreatic) VET test

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







English

ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C)!

1. Preparing the test system

- 1.1 Insert RFID card
- 1.2 Place ERS cuvette into test kit rack
- 1.3 Place ERS cap into test kit rack
- 1.4 Press "Measurement" button and enter all required information using the touchscreen

2. Sample preparation

Choose sample type:

- 2.1a "Pipette20µ", measurement range: 25 - 300 U/L
OR ...
- 2.1b "Pipette5µ", measurement range: 250 - 1200 U/L
- 2.2 Aspirate 20 µl **OR** 5 µl sample material from centrifuged sample tube

3. Test processing

- 3.1 Dispense 20 µl **OR** 5 µl sample INTO THE LIQUID in the ERS cuvette
- 3.2 Firmly apply ERS cap on ERS cuvette
- 3.3 Place ERS cartridge into the solo analyser
- 3.4 Start automatic sample processing by pressing the  start button on the solo analyser