# NH<sub>3</sub> Ammonia VET test kit

For veterinary use only!

Veterinary test kit for quantitative in vitro determination of NH<sub>3</sub> Ammonia in EDTA plasma with a solo analyser.

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**Test kit preparation: allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Return the test kit package to the refrigerator.**

## Summary

Protein digestion produces ammonia within the intestine by deamination of amino acid. Ammonia reaches the liver by portal blood where it is synthesized into urea (main route of metabolism of intestinal ammonia into urea).

Heavy liver dysfunction (e.g. shunt, cirrhosis) can prevent the liver from synthesizing ammonia into urea. As a result one will receive a lowered urea concentration in relation to creatinine as well as increased blood ammonia. Ammonia and other protein metabolism (aromatic amino acids) that get out of control of the liver metabolites can affect the central nervous system (hepatoencephalic syndrome).

## Method

Enzymatic test. Photometric measurement of the absorbance rate (kinetic) at 340 nm wavelength.

## Measurement Range

10 - 1100 µmol/l (Lot dependent)

## Sample Material

Use only 20 µl fresh EDTA plasma (Ammonia heparin plasma or serum must NOT be used)! After drawing the blood, immediately centrifuge and analyse within 5 minutes on the analyser.

## Test Kit

ERS cuvette filled with buffer, pH 8.0, alpha-ketoglutarate, GLDH and NADPH.

## Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 - 8 °C.  
Opening the vacuum package may limit the reagent stability to three months (stored at 2 - 8 °C) from the date of opening.  
DO NOT FREEZE!

## Warnings and Precautions

DO NOT INGEST! Avoid contact with skin and eyes. Observe all necessary precautions for the use of laboratory reagents.

## Waste Management

Always refer to local legal requirements.

## Reference Range

Dog: < 60 µmol/l  
Cat: < 60 µmol/l

It is recommended that each laboratory establishes its own reference ranges.

## Quality Control

For internal quality control the Eurolyser NH<sub>3</sub> Ammonia VET control kit is recommended. Order number: VT 1600

## Precision

Reproducibility within-run:  
Control; N = 20; mean = 246 µmol/l; CV = 4.9%;

## Correlation

$y$  (solo NH<sub>3</sub>) = 0.972 x (reference method NH<sub>3</sub>) + 0.233; R<sup>2</sup> = 0.985;

## Interferences

The test system has been analysed for various interferences. Criterion was the recovery within 10% of initial values.

Ascorbic acid	3 mg/dl
Bilirubin	40 mg/dl
Lactate	200 mg/dl
Pyruvate	6.6 mg/l
Triglycerides	200 mg/dl

## References

- HARDY, RM. (1989) Diseases of the liver and their treatment, In: EttingerSJ, ed. Textbook of Veterinary Internal Medicine, 3rd ed. Philadelphia : WB Saunders; 1479 – 1527.
- RASKIN, NH., PRICE, JB., FISHMAN, RA. (1964) Portal-systemic encephalopathy due to congenital intrahepatic shunts. N Engl J Med; 270:225 – 229;
- RASKIN, NH., BREDESEN, D., EHRENFIELD, WK., et al. (1984) Periodic confusion caused by congenital extrahepatic portacaval shunt; Neurology; 34:666 – 669;

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
## References


- HARDY, RM. (1989) Diseases of the liver and their treatment, In: EttingerSJ, ed. Textbook of Veterinary Internal Medicine, 3rd ed. Philadelphia : WB Saunders; 1479 – 1527.
- RASKIN, NH., PRICE, JB., FISHMAN, RA. (1964) Portal-systemic encephalopathy due to congenital intrahepatic shunts. N Engl J Med; 270:225 – 229;
- RASKIN, NH., BREDESEN, D., EHRENFELD, WK., et al. (1984) Periodic confusion caused by congenital extrahepatic portacaval shunt; Neurology; 34:666 – 669;

# Durchführung eines NH<sub>3</sub> Ammoniak VET Tests Performing of a NH<sub>3</sub> Ammonia VET test

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





---

### **ATTENTION!**

Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C)!

### **1. Preparing the test system**

- 1.1 Insert RFID card
- 1.2 Place ERS cuvette into test kit rack
- 1.3 Detach cap from ERS cuvette
- 1.4 Press "Measurement" button and enter all required information using the touchscreen

### **2. Sample preparation**

- 2.1 Aspirate 20 µl sample material from centrifuged sample tube

### **3. Test processing**

- 3.1 Dispense 20 µl sample INTO THE LIQUID in the ERS cuvette
- 3.2 Firmly apply ERS cap on ERS cuvette
- 3.3 Place ERS cartridge into the solo analyser
- 3.4 Start automatic sample processing by pressing the  start button on the solo analyser