Phenobarbital VET test kit

For veterinary use only!

Veterinary test kit for quantitative in vitro determination of phenobarbital in serum or lithium heparin plasma with a solo analyser.

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Order number: VT0211
Order number: VT2100

Indication

Phenobarbital VET test kit
Phenobarbital VET control kit

Kit size

6 tests/kit
1 x 5 ml (decision level)

Test kit preparation: allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Return the test kit package to the refrigerator.

Summary

Phenobarbital is administered to control seizures and epilepsy in dogs and cats. As a barbiturate it loosens cramps. Depending on the dosage, a sedative or hypnotic action of the drug can be observed. In cats, phenobarbital is the first choice as anti-epileptic treatment. In dogs, it may be administered together with potassium-bromide, depending on individual reaction of the patient to the treatment.

Blood phenobarbital concentration is recommended to be checked approximately 12 - 14 days after administration. Best sampling time is approximately two hours prior to the next scheduled drug treatment.

Method

Quantitative turbidimetric test for the measurement of phenobarbital in serum and lithium-heparin plasma. Latex particles coated with phenobarbital are agglutinated when mixed with phenobarbital antibody solution. When a sample containing phenobarbital is used, the agglutination reaction is partially inhibited, slowing the process of agglutination. The agglutination causes an absorbance change (546 nm), inversely dependent upon the phenobarbital contents of the patient sample.

Measurement Range

Measurement range: 0 - 100 µg/ml (Lot dependent)
Sensitivity: 10 µg/ml

Sample Material

20 µl serum (preferred) or lithium-heparin plasma.

Sample material stability:

Stored at 2 - 8 °C: 3 days
Stored at - 20 °C: 60 days

Test Kit

ERS cuvette filled with buffer reagent
ERS cap filled with phenobarbital sensitized latex

Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 - 8 °C.

Opening the vacuum package may limit the reagent stability to three months (stored at 2 - 8 °C) from the date of opening. DO NOT FREEZE!

Warnings and Precautions

DO NOT INGEST! Avoid contact with skin and eyes. Observe all necessary precautions for the use of laboratory reagents.

Waste Management

Always refer to local legal requirements.

Reference Range

Dog: 15 - 35 µg/ml¹
Cat: 15 - 35 µg/ml²
Horse: 15 - 45 µg/ml³
Foal: 5 - 30 µg/ml

It is recommended that each laboratory establishes its own reference ranges.

Quality Control

For internal quality control the Eurolyser Phenobarbital VET control kit is recommended. Order number: VT 2100

Precision

Reproducibility within-run:
Control Serum; N = 12; mean = 20.0 µg/ml; CV = 9.98%
Control Serum; N = 12; mean = 45.1 µg/ml; CV = 6.10%

Correlation

30 samples splitted in serum and lithium-heparin:
 y (solo PHENOB) = $0.8874 x$ (Immulite 2000 PHENOB) + 1.1693;
 $R^2 = 0.8502$;

Interferences

The following substances show no interferences up to:

Bilirubin	25 mg/dl
Triglycerides	1000 mg/dl
Haemoglobin	1000 mg/dl
Intralipid	800 mg/dl

References

¹<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4913615/>

RAVIS WR., PEDERSOLI WM., WIKE JS. (1989) Pharmacokinetics of phenobarbital in dogs given multiple doses. Am J Vet Res;50:1343-1347. [PubMed]

SCHWARTZ-PORSCHÉ D., LOSCHER W., FREY HH. (1985) Therapeutic efficacy of phenobarbital and primidone in canine epilepsy: a comparison. J Vet Pharmacol Therap;8:113-119. [PubMed]

² THOMAS WB., DEWEY CW. (2003) Seizures and narcolepsy: A Practical Guide to Canine and Feline Neurology—Ames: Iowa State University Press, pp 193-212.

BOOTHE DM. (2000) Anticonvulsants and other neurologic therapies in small animals: Small Animal Clinical Pharmacology and Therapeutics, 1st ed— Philadelphia: WB Saunders, pp 431-456.

³ FURR M., REED JARVIS C S. (1976) Equine Neurology, The evolutionary strategy of the equidae and the origins of rumen and cecal digestion. Evolution 30:757-74.

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SCHWARTZ-PORSCHKE D., LOSCHER W., FREY HH. (1985) Therapeutic efficacy of phenobarbital and primidone in canine epilepsy: a comparison. J Vet Pharmacol Therap;8:113-119. [PubMed]

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
Performing a **Phenobarbital** VET test

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





English

ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C)!

1. Preparing the test system

- 1.1 Insert RFID card
- 1.2 Place ERS cuvette into test kit rack
- 1.3 Place ERS cap into test kit rack
- 1.4 Press "Measurement" button and enter all required information using the touchscreen

2. Sample preparation

- 2.1 Aspirate 20 µl sample material from centrifuged sample tube

3. Test processing

- 3.1 Dispense 20 µl sample INTO THE LIQUID in the ERS cuvette
- 3.2 Firmly apply ERS cap on ERS cuvette
- 3.3 Place ERS cartridge into the solo analyser
- 3.4 Start automatic sample processing by pressing the  start button on the solo analyser