SAA VET test kit

For veterinary use only!

Veterinary test kit for quantitative in vitro determination of SAA in serum with a solo analyser.

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Kit size

6 tests/kit
1 x 2 ml

Test kit preparation: allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Return the test kit package back to the refrigerator.

Summary

Serum amyloid A (SAA) is an early and sensitive blood biomarker for tissue injury and inflammation and has been indicated in many inflammatory diseases. The level of SAA circulating in the blood is known to increase dramatically in response to tissue damage or inflammation, classifying it as an acute phase protein. Circulating SAA concentrations may increase up to 1000-fold following inflammation, infection, tissue injury and cell necrosis and decline rapidly following recovery.

It has also been observed in several veterinary species including the horse that SAA is a very useful inflammatory marker that may be used for detection of clinical and possibly subclinical disease, monitoring of disease activity and response to therapy, prognostication, and detection of spread of infections in herds (Nielsen et al., 2004). The kinetic profile of the SAA response makes SAA an excellent indicator of inflammation. Hepatic SAA synthesis begins shortly after an inflammatory insult. As a result plasma concentrations start to increase within few hours, and they peak 36 – 48 h after injury, as demonstrated previously in horses after experimental induction of inflammatory or infectious disease (Hulte'n et al., 2002).

Method

Photometric endpoint measurement at 546 nm wavelength.

Latex agglutination test based on optical measurement of the change in turbidity caused by the agglutination of the latex particles sensitised with SAA antibodies.

Measurement Range

Horse: 10 - 500 µg/ml (10 - 500 mg/l) (Lot dependent)
Cat: 10 - 150 µg/ml (10 - 150 mg/l) (Lot dependent)

Samples with concentrations higher than the upper limit of the measurement range must be diluted 1 + 5 with physiological saline (0.9% NaCl), e.g. 20 µl sample + 100 µl 0.9% NaCl solution, and the result then multiplied by 6.

Sample Material

Use serum or alternatively lithium-heparin plasma.

Test Kit

ERS cuvette filled with Good's buffer.

ERS cap filled with anti-human SAA rabbit polyclonal antibodies and anti-human SAA mouse monoclonal antibodies.

Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 - 8 °C.

Opening the vacuum package may limit the reagent stability to three months (stored at 2 - 8 °C) from the date of opening.

DO NOT FREEZE!

Warnings and Precautions

DO NOT INGEST! Avoid contact with skin and eyes. Contains sodium azide, which may react with lead or copper plumbing to form explosive compounds. Observe all necessary precautions for the use of laboratory reagents.

Waste Management

Always refer to local legal requirements.

Reference Range

Cat: < 10 µg/ml
Horse: < 20 µg/ml
Camel: < 10 µg/ml

It is recommended that each laboratory establishes its own reference ranges.

Quality Control

For internal quality control the Eurolyser SAA VET control kit is recommended. Order number: VT 1700

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English

ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C)!

1. Preparing the test system

- 1.1 Insert RFID card
 - 1.2 Place ERS cuvette into test kit rack
 - 1.3 Place ERS cap into test kit rack
 - 1.4 Press "Measurement" button and enter all required information using the touchscreen
-

2. Sample preparation

- 2.1 Aspirate 5 µl sample material from centrifuged sample tube
-

3. Test processing

- 3.1 Dispense 5 µl sample INTO THE LIQUID in the ERS cuvette
 - 3.2 Firmly apply ERS cap on ERS cuvette
 - 3.3 Place ERS cartridge into the solo analyser
 - 3.4 Start automatic sample processing by pressing the  start button on the solo analyser.
-