T4-VET test kit

For veterinary use only!

Veterinary test kit for quantitative in vitro determination of T4 with solo analyser.

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Kit size

16 tests/pack
2 x 1 ml (low/high)

Test kit preparation: Allow single test min. 10 minutes to warm up to room temperature by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

Thyroxine (T4) is produced by the thyroid. Controlled by a complex regulatory mechanism it is released into the bloodstream. Most of the thyroxine in blood is bound to various carrier proteins such as albumine, prealbumine and other specific binding proteins. Diseases of the thyroid with abnormal hormone concentrations occur more likely in older animals. Elevated or lowered concentrations of the thyroid hormones, e.g. thyroxine (T4), lead towards problems in numerous other systems. Hyperthyreosis are merely found in cats, hypothyreosis are found in all species.

Method

The T4-VET test is a homogenous enzymatic immunoassay. By using an 8-anilino-1-naphtaline acid sulfure (ANS) the bound T4 is seperated from its binding proteins. The seperated thyroxine released from blood competes afterwards with a thyroxine that is linked to the Glucose-6-phosphat dehydrogenase enzyme to a specified binding area. The GAPDH activity at 340 nm is associated with the complete thyroxine concentration in blood. If no T4 is in the blood sample the enzyme activity is inhibited by the binding of a specific antibody.

Measurement Range

When using a 20 µl pipette:

0 - 8 µg/dl (0 - 103 nmol/l)

When using a 5 µl pipette:

8 - 32 µg/dl (103 - 413 nmol/l)

Sample material

Serum or plasma (EDTA, Heparin, Citrat).

Samplevolume according to the menu settings on the laboratory photometer and according to the use of a 20 µl or 5 µl pipette.

Stability and Storage

Store at 2 - 8°C. DO NOT FREEZE. Kit can be used until expiry date printed on the labels.

Warnings

Do NOT INGEST. Avoid contact with skin and eyes. Contains sodium azide, which may react with lead or copper plumbing to form explosive compounds. Take the necessary precautions for the use of laboratory reagents.

Normal Ranges

Dog: 1,3 - 4,5 µg/dl (17 - 58 nmol/l)

Cat: 1,0 - 4,0 µg/dl (13 - 52 nmol/l)

Pets: 0,8 - 2,5 µg/dl (10 - 32 nmol/l)

Horse: 1,3 - 4,1 µg/dl (17 - 53 nmol/l)

Cattle: 3,8 - 8,2 µg/dl (49 - 106 nmol/l)

Precision

Sample 20 µl:

Dog; rep.=10; mean=1,03µg/dl; SD=0,04µg/dl; CV=3,8%

Cat; rep.=10; mean=1,93µg/dl; SD=0,08µg/dl; CV=4,1%

Cat; rep.=14; mean=7,98µg/dl; SD=0,35µg/dl; CV=4,4%

Sample 5 µl:

Cat; rep.=5; mean=7,79µg/dl; SD=0,27µg/dl; CV=3,5%

Correlation

Method comparison with Immulite 1000/canine Total T4:
Eurolyser solo T4-VET (µg/dl) = 0,978x + 0,116; R² = 0,965;
N = 43 cats;

Method comparison with Immulite 1000/canine Total T4:
Eurolyser solo T4-VET (µg/dl) = 0,990x - 0,032; R² = 0,951;
N = 48 dogs;

Method comparison with Immulite 1000/canine Total T4:
Eurolyser solo T4-VET (µg/dl) = 0,985x + 0,035; R² = 0,963;
N = 48 dogs and 43 cats; mean = 2,38 µg/dl

Waste Management

Please refer to local legal requirements.

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English

WARNING!

Do not wipe pipetting tip after sample collection! Allow test minimum 10 minutes to warm up to room temperature.

1. Prepare testsystem

- 1.1 Insert RFID card
- 1.2 Place ERS cuvette into testkit rack.
- 1.3 Place ERS cap into testkit rack.
- 1.4 Press measurement button, enter required information using the touch screen

2. For further details please see to user manual of laboratory photometer

Choose sample type:

- 2.1a 20 µl pipette: Possible measuring range:
0 - 8 µg/dl (0 - 103 nmol/l)
OR...
- 2.1b 5 µl pipette: Possible measuring range:
8 - 32 µg/dl (103 - 413 nmol/l)
- 2.2 Start of analysis by pressing the start button

3. Sample preparation serum or plasma (serum primary tube, lithium-heparin or EDTA plasma)

- 3.1 Pipette 20 µl or 5 µl serum/plasma from primary tube
- 3.2 Pipette 20 µl or 5 µl sample into ERS cuvette INTO LIQUID
- 3.3 Place ERS cap firmly on ERS cuvette
- 3.4 Place ERS cartridge into laboratory photometer
- 3.5 Start automatic sample processing by pressing the  start button.