Phenobarbital VET test kit **EURO**Lyser

For veterinary use only!

Veterinary test kit for quantitative in vitro determination of phenobarbital in serum or lithium heparin plasma on a solo or CUBE-VET analyser

 Eurolyser Diagnostica GmbH
Bindergasse 3
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00
Fax +43 (0)662 / 43 21 00-50
www.eurolyser.com

English

Order information

Order number: VT 0211
Order number: VT 2100

Indication

Phenobarbital VET test kit
Phenobarbital VET control kit

Kit size

6 tests/pack
1 x 5 ml (decision level)

Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 – 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

Phenobarbital is administered to control seizures and epilepsy in dogs and cats. As a barbiturate it loosens cramps. Depending on the dosage, a sedative or hypnotic action of the drug can be observed. In cats, phenobarbital is the first choice as anti-epileptic treatment. In dogs, it may be administered together with potassium bromide, depending on individual reaction of the patient to the treatment.

Blood phenobarbital concentration is recommended to be checked approximately 12 – 14 days after administration. The best sampling time is about two hours prior to the next scheduled drug treatment.

Method

Homogeneous, immunoturbidimetric test.

Measurement Range

5 – 80 µg/ml

Sample Material

Use 20 µl of serum or lithium-heparin plasma.

Test Kit

ERS cuvette filled with buffer reagent
ERS cap filled with antibody reagent

Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 – 8 °C. Opening the vacuum package may limit the reagent stability to three months (stored at 2 – 8 °C) from the date of opening. **DO NOT FREEZE!**

Warnings and Precautions

DO NOT INGEST! Avoid contact with skin and eyes. Observe all necessary precautions for the use of laboratory reagents.

Waste Management

Please refer to local legal requirements.

Reference Range

Dog: 15 – 35 µg/ml [1]
Cat: 15 – 35 µg/ml [2]
Horse: 15 – 45 µg/ml [3]
Foal: 5 – 30 µg/ml

It is recommended that each laboratory establishes its own reference ranges.

Quality Control

For internal quality control the Eurolyser Phenobarbital VET control kit is recommended. Order number: VT2100

Precision

Reproducibility within-run:
Control; N = 20; mean = 26.7 µg/ml; CV = 7.62%;

Correlation

Sample correlation: N = 75

$y(\text{Eurolyser Phenob}) = 0.9563x(\text{Immulite Phenob}) + 1.1396; R^2 = 0.9269$

Interferences

The test system has been analysed for various interferences. Criterion was the recovery within 10% of initial values.

Haemoglobin	525 mg/dl
Human albumin	12 g/dl
Bilirubin conjugated	36 mg/dl
Bilirubin unconjugated	36 mg/dl
Cholesterol	620 mg/dl
Rheumatoid factor	1080 IU/ml
Triglycerides	835 mg/dl
Uric acid	30 mg/dl

References

- [1] RAVIS WR., PEDERSOLI WM., WIKE JS. (1989) Pharmacokinetics of phenobarbital in dogs given multiple doses; *Am J Vet Res*; 50:1343 – 1347;
SCHWARTZ-PORSCHKE D., LOSCHER W., FREY HH. (1985) Therapeutic efficacy of phenobarbital and primidone in canine epilepsy: a comparison; *J Vet Pharmacol Therap*; 8:113 – 119;
- [2] THOMAS WB., DEWEY CW. (2003) Seizures and narcolepsy: A practical guide to canine and feline neurology - Ames: Iowa state university press, pp. 193 – 212;
BOOTHE DM. (2000) Anticonvulsants and other neurologic therapies in small animals: *Small Animal Clinical Pharmacology and Therapeutics*, 1st ed - Philadelphia: WB Saunders, pp. 431 – 456;
- [3] FURR M., REED JARVIS CS. (1976) Equine neurology, the evolutionary strategy of the equidae and the origins of rumen and cecal digestion. *Evolution* 30:757 – 74;

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Referenzen

- [1] RAVIS WR., PEDERSOLI WM., WIKE JS. (1989) Pharmacokinetics of phenobarbital in dogs given multiple doses; Am J Vet Res; 50:1343 – 1347;
SCHWARTZ-PORSCHKE D., LOSCHER W., FREY HH. (1985) Therapeutic efficacy of phenobarbital and primidone in canine epilepsy: a comparison; J Vet Pharmacol Therap; 8:113 – 119;
- [2] THOMAS WB., DEWEY CW. (2003) Seizures and narcolepsy: A practical guide to canine and feline neurology - Ames: Iowa state university press, pp. 193 – 212;
BOOTHE DM. (2000) Anticonvulsants and other neurologic therapies in small animals: Small Animal Clinical Pharmacology and Therapeutics, 1st ed - Philadelphia: WB Saunders, pp. 431 – 456;
- [3] FURR M., REED JARVIS CS. (1976) Equine neurology, the evolutionary strategy of the equidae and the origins of rumen and cecal digestion. Evolution 30:757 – 74;

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Processing of a **Phenobarbital** VET test

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English

ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 – 25 °C) before use!

1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Place ERS cuvette in test kit rack
- 1.3 Place ERS cap in test kit rack
- 1.4 Press „Measurement“ button, enter required information using the touchscreen

2. Sample preparation

- 2.1 Aspirate 20 µl sample material from centrifuged sample tube

3. Sample processing

- 3.1 Dispense 20 µl sample INTO THE LIQUID in the ERS cuvette
- 3.2 Apply ERS cap firmly onto ERS cuvette
- 3.3 Place ERS cartridge into laboratory photometer
- 3.4 Start automatic sample processing by pressing the  start button on the solo laboratory photometer, or by closing the door of the CUBE-VET laboratory photometer.