# Lactate VET test kit

For veterinary use only!

Veterinary test kit for quantitative in vitro determination of lactate in plasma with a solo analyser.

**EURO**Lyser 

 Eurolyser Diagnostica GmbH Tel. +43 (0)662 / 43 21 00  
Bayernstraße 11a Fax +43 (0)662 / 43 21 00-50  
5020 Salzburg/Austria www.eurolyser.com

## English

### Order information

Order number: VT 0220  
Order number: VT 0221  
Order number: VT 2200

### Indication

Lactate VET test kit  
Lactate VET test kit  
Lactate VET control kit

### Kit size

16 tests/pack  
6 tests/pack  
1 x 5 ml (decision level)

**Allow single test to warm up at room temperature (20 - 25 °C) by placing the test into the test kit rack for a minimum of 10 minutes before use. Put the test kit package back into refrigerator.**

### Summary

In the face of anaerobic metabolism lactate is synthesized from pyruvate. Muscle tissue is the predominant source of lactate production. Increased lactate concentration is a sensitive indication for a significant metabolic disturbance. Possible reasons are hypoxia and hypoperfusion. Diseases associated with increased lactate concentration are i.e. hypovolemic shock, local hypoperfusion (i.e. intestinal ischaemia in the face of ileus, aortic-thromboembolism), hypoxia, anaemia, or sepsis. In racehorses lactate is used to either check performance ability or i.e. exclude latent myopathy.

### Method

Colorimetric test. Photometric endpoint measurement of the absorbance at 546 nm wavelength.

### Measurement Range

0.1 – 15 mmol/l

When values exceed this range, samples should be diluted 1 + 1 with NaCl solution (9 g/l) and the results multiplied by 2.

### Sample Material

Use 20 µl of plasma. Do not use serum!

For measurements right after taking the blood sample:  
Use li-heparin samples and separate plasma from red blood cells immediately (within 5 minutes).

In case measurements cannot be done immediately:  
Use Na-fluorid/EDTA samples, centrifuge and measure at the earliest.

Stability in plasma (after centrifugation and separation):  
8 hours at 15 – 25 °C, 14 days at 2 – 8 °C.

### Test Kit

ERS cuvette filled with mono reagent.

### Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 – 8 °C.  
Opening the vacuum package may limit the reagent stability to 3 months (stored at 2 – 8 °C) from the date of opening. DO NOT FREEZE!

### Warnings and Precautions

DO NOT INGEST! Avoid contact with skin and eyes. Observe all necessary precautions for the use of laboratory reagents.

### Waste Management

Always refer to local legal requirements.

### Reference Range

Canine: 0.5 – 3.0 mmol/l  
Feline: < 1.0 mmol/l  
Equine: 0.5 – 2.0 mmol/l

It is recommended that each laboratory establishes its own reference ranges.

### Quality Control

For internal quality control the Eurolyser Lactate VET control kit is recommended. Order number: VT 2200

### Precision

Reproducibility within-run:

Control high; N = 20; mean = 3.87 mmol/l; CV = 1.69%;  
Control low; N = 20; mean = 1.49 mmol/l; CV = 2.10%;

### Correlation

Equine sample correlation: N = 40

$y$  (solo Lactate) = 1.0602 x (Cobas Lactate) - 0.2564;  
 $R^2$  = 0.9577;

### Interferences

No interference up to:

Ascorbic acid	30 mg/dl
Bilirubin	60 mg/dl
Triglyceride	2000 mg/dl
Haemoglobin	1000 mg/dl
Dopamine	10 mg/l
L-Dopamine	20 mg/l
Methyldopamine	10 mg/l
Glycolic acid	1200 mg/l

For further information on interfering substances refer to Young DS<sup>1</sup>.

### References

1. YOUNG, D.S. (2000) Effects of Drugs on Clinical Laboratory Tests. 5th ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press;

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Bestellnummer: VT 0220  
Bestellnummer: VT 0221  
Bestellnummer: VT 2200

### Indication

Lactate VET test kit  
Lactate VET test kit  
Lactate VET control kit

### Kit size

16 Tests/Packung  
6 Tests/Packung  
1 x 5 ml (decision level)

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---

### ATTENTION!

Allow single test to warm up to room temperature (20 - 25 °C) for a minimum of 10 minutes before use!

- 1. Preparation of test system**
  - 1.1 Insert RFID card
  - 1.2 Place ERS cuvette into test kit rack
  - 1.3 Take off ERS cuvette cap
  - 1.4 Press „Measurement“ button and enter all required information using the touchscreen

- 
- 2. Sample preparation**
    - 2.1 Aspirate 20 µl sample material from centrifuged sample tube

- 
- 3. Test processing**
    - 3.1 Dispense 20 µl sample INTO THE LIQUID in the ERS cuvette
    - 3.2 Firmly apply cap onto ERS cuvette
    - 3.3 Place ERS cartridge into the solo analyser
    - 3.4 Start automatic sample processing by pressing the  start button on the solo analyser