# Cortisol

## VET test kit

For veterinary use only!

Intended use:

Veterinary test kit for quantitative in vitro determination of cortisol in serum on a solo or CUBE-VET analyser

EUROlyser Diagnostica GmbH  
Bindergasse 3  
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00  
Fax +43 (0)662 / 43 21 00-50  
www.eurolyser.com

## English

### Order information

Order number: VT0290  
Order number: VT0291  
Order number: VT2900

### Indication

Cortisol VET test kit  
Cortisol VET test kit  
Cortisol VET control kit

### Kit size

16 tests  
6 tests  
1 x 0.5 ml (decision level)

**Important: AFTER preparing the reagent, allow single test at least 10 min. to warm up to room temperature (20 – 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.**

### Summary

Cortisol is a glucocorticoid hormone that affects the metabolism of carbohydrates, proteins, and fats. Chronic, excessive glucocorticoid release leads to hyperadrenocorticism (HAC). Pituitary-dependent hyperadrenocorticism (PDH), more commonly known as Cushing's disease, is caused by a pituitary tumour (mostly adenomas) in the brain. The tumour triggers excessive levels of adrenocorticotropic hormone (ACTH), which stimulates the adrenal glands near the kidneys to produce cortisol. A small percentage of dogs with HAC have a tumour of one of the adrenal glands. This form of HAC is called adrenal-dependent HAC and results from a direct increase in cortisol production by the adrenal gland tumour. The overproduction of cortisol causes symptoms such as hair loss, pot-bellied appearance, increased appetite, and polydipsia and polyuria. Addison's disease (hypoadrenocorticism, or underproduction of cortisol) occurs less commonly than the Cushing's disease in dogs.

### Method

Homogeneous, immunoturbidimetric test

### Measurement Range

0.8 – 24 µg/dl

### Sample Material

Use 40 µl of serum

### Test Kit

ERS cuvette filled with buffer reagent  
ERS cap filled with antibody reagent  
Bottle filled with dilution buffer

### Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 – 8 °C. Opening the vacuum package may limit the reagent stability to 3 months (stored at 2 – 8 °C) from the date of opening. DO NOT FREEZE!

### Warnings and Precautions

DO NOT INGEST! Avoid contact with skin and eyes. Observe all necessary precautions for the use of laboratory reagents.

### Waste Management

Please refer to local legal requirements.

### Interferences

The test system has been analysed for various interferences. Criterion was the recovery within 15% of initial values.

Haemoglobin	131 mg/dl
Albumin	12 g/dl
Bilirubin, conjugated	36 mg/dl
Bilirubin, unconjugated	36 mg/dl
Cholesterol	620 mg/dl
Rheumatoid factor	540 IU/ml
Triglycerides	835 mg/dl
Uric acid	30 mg/dl

### Precision

Reproducibility within-run:  
Control; N = 20; mean = 5.7 µg/dl; CV = 6.93%

### Correlation

Sample correlation: N = 47

y (Eurolyser cortisol) = 0.9168x (Roche Elecsys) - 0.2049; R<sup>2</sup> = 0.9624;

### Quality Control

For internal quality control the Eurolyser Cortisol VET control kit is recommended. Order number: VT2900

### Reference Ranges

It is recommended that each laboratory establishes its own reference ranges.

#### I. ACTH Stimulation Test

Pre-ACTH	Post-ACTH	Interpretation
< 2 µg/dl	< 2 µg/dl	Consistent with hypoadrenocorticism
	2 – 6 µg/dl	Inconclusive
2 – 6 µg/dl	6 – 18 µg/dl	Normal
	18 – 22 µg/dl	Equivocal, Cushing's syndrome possible
> 22 µg/dl	> 22 µg/dl	Consistent with Cushing's syndrome

#### II. Low-dose Dexamethasone Suppression Test

4-hour cortisol level	8-hour cortisol level	Interpretation
< 1 µg/dl	< 1 µg/dl	Normal
1 – 1.5 µg/dl	1 – 1.5 µg/dl	Inconclusive
> 1.5 µg/dl and > 50% of baseline	> 1.5 µg/dl and > 50% of baseline	Consistent with Cushing's syndrome
< 1.5 µg/dl or < 50% of baseline	> 1.5 µg/dl and > 50% of baseline	Consistent with PDH
> 1.5 µg/dl or > 50% of baseline	> 1.5 µg/dl and < 50% of baseline	Consistent with PDH

#### III. High-dose Dexamethasone Suppression Test

4-hour cortisol level	8-hour cortisol level	Interpretation
< 1.5 µg/dl or < 50% of baseline	> 1.5 µg/dl and > 50% of baseline	Consistent with PDH
> 1.5 µg/dl and > 50% of baseline	< 1.5 µg/dl or < 50% of baseline	Consistent with PDH
< 1.5 µg/dl or < 50% of baseline	< 1.5 µg/dl or < 50% of baseline	Consistent with PDH
> 1.5 µg/dl and > 50% of baseline	> 1.5 µg/dl and > 50% of baseline	Consistent with PDH or ATH

#### IV. Post ACTH Stimulation Mitotane (Lysodren®) Dosing and Monitoring

Loading Dose	Interpretation	Maintenance Dose	Interpretation
< 1 µg/dl	Discontinue mitotane	< 1 µg/dl	Discontinue mitotane
1 – 5 µg/dl	Begin maintenance mitotane dosing	1 – 5 µg/dl	Maintain current dosage
> 5 µg/dl	Continue mitotane loading dose	> 5 µg/dl	Increase weekly dose or repeat loading dose.

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## English

### ATTENTION!

**AFTER** preparing the reagent, allow single test at least 10 minutes to warm up to room temperature (20 – 25 °C) before use!

#### 1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Press „Measurement“ button, enter required information using the touchscreen
- 1.3 Place ERS cuvette in test kit rack
- 1.4 Place ERS cap in test kit rack
- 1.5 Take out bottle containing dilution buffer from the test kit package

#### 2. Reagent and sample preparation

- 2.1 Aspirate 400 µl from the dilution buffer bottle
- 2.2 Dispense 400 µl dilution buffer INTO THE LIQUID in the ERS cuvette

**Allow single test to warm up to room temperature for at least 10 minutes before measurement! Put test kit package back into refrigerator**

- 2.3 Aspirate 40 µl sample material from centrifuged sample tube

#### 3. Sample processing

- 3.1 Dispense 40 µl sample INTO THE LIQUID in the ERS cuvette
- 3.2 Apply ERS cap firmly onto ERS cuvette
- 3.3 Place ERS cartridge into laboratory photometer
- 3.4 Start automatic sample processing by pressing the  start button on the solo laboratory photometer, or by closing the door of the CUBE-VET laboratory photometer.