

ACUTE CARE TESTING

ABL90 FLEX analyzer



Time on your side

Now time is on *your* side

ABL90 FLEX 17 parameters in 35

Fast results and always ready

The ABL90 FLEX analyzer offers automatic sample mixing in just seven seconds and full acute care panel results in only 35 seconds.

Easy to use

Reduce the risk of user errors and save time.

Easy to maintain

Eliminate maintenance surprises and ensure high uptime.
Easy steps to be performed by any user at any time of the day or night.

ABL90 FLEX analyzer sensor



Test menu

pH | $p\text{CO}_2$ | $p\text{O}_2$ | $s\text{O}_2$ | ctHb | FO_2Hb | FCOHb | FMetHb | FHHb
 FHbF | cK^+ | cNa^+ | cCa^{2+} | cCl^- | cGlu | cLac | ctBil



EX analyzer

5 seconds on 65 μ L



Meets regulatory compliance

Compliance to regulatory requirements with complete data integrity

Accuracy and pre-analytical error reduction

Get reliable, accurate and precise results by reducing patient/sample mix-up and securing consistent sample mixing

1st Automatic ready

Get the right result for the right patient at the right time, while improving patient and operator safety and sample integrity. 1st Automatic simplifies your blood gas testing from test ordered to reported results.

