



Samsung PT10V - Testkassetter

Parameter	Electrolyte Test	Comprehensive Test	Diabetes Test	Pre surgical Test	Kidney Test	Kidney Monitoring Test	Liver Test	Liver Monitoring Test	Large Animal Test	Parameter
A/G ratio*		X		X			X		X	A/G ratio*
ALB		X		X	X		X		X	ALB
ALP		X		X			X		X	ALP
ALT/GPT		X		X			X			ALT/GPT
AMYL		X			X					AMYL
AST/GOT							X		X	AST/GOT
B/C ratio*		X		X	X	X			X	B/C ratio*
BUN		X		X	X	X	X		X	BUN
CA		X			X					CA
CHOL		X	X							CHOL
CK									X	CK
CREA		X		X	X	X			X	CREA
Fruktosamin			X							Fruktosamin
GGT							X		X	GGT
GLOB*		X		X			X		X	GLOB*
GLU		X	X	X	X		X		X	GLU
PHOS		X			X	X				PHOS
NA <sup>+</sup> / K <sup>+</sup> / Cl <sup>-</sup>	X									NA <sup>+</sup> / K <sup>+</sup> / Cl <sup>-</sup>
Na/K ratio*	X									Na/K ratio*
TBIL		X					X		X	TBIL
TG		X	X							TG
TP		X		X			X		X	TP

2018-06-18



\* Beräknad parameter



\* Beräknad parameter