

Fuji NX500i - paneler och singeltester											
Parameter		Singeltester		Comprehensive panel		Pre surgery panel		Kidney Profile panel		Liver panel	Parameter
ALB		X		X				X		X	ALB
ALP		X		X		X				X	ALP
ALT/GPT		X		X		X				X	ALT/GPT
AMYL		X									AMYL
AST/GOT		X									AST/GOT
BUN		X		X		X		X			BUN
CA		X		X				X			CA
CHOL		X		X							CHOL
CPK		X									CK
CREA		X		X		X		X			CREA
cCRP		X									cCRP
GGT		X		X						X	GGT
GLU		X		X		X				X	GLU
IP/PHOS		X		X				X			IP/PHOS
NA ⁺ /K ⁺ /Cl ⁻		X									NA ⁺ /K ⁺ /Cl ⁻
LDH		X									LDH
LIPAS		X									LIP
MG ⁺⁺		X									MG ⁺⁺
NH ₃		X									NH ₃
TBIL		X		X						X	TBIL
TG		X									TG
TP		X		X		X		X			TP
UA		X									UA
Test/fp		24		6x12		12x6		4x6		4x6	Test/fp



 TRIOLAB 

Tel: 031-81 72 00
www.triolab.se
info@triolab.se