

Rotorer till Vetscan Classic och VS2

Rotor	Avian / Reptilian Plus*	Canine Wellness	Comprehensive Diagnostic	Critical Care Plus*	Electrolyte Plus -- NYHET!	Equine Plus*	Kidney Profile Plus	Large Animal*	Mammalian Liver*	Prep Profile II*	Preventive Care Profile Plus -- NYHET!	Thyroxine / Cholesterol*	Phenobarbital Profile -- NYHET!	Parameter
ALB	x	x	x			x	x	x	x		x			ALB
ALP		x	x					x	x	x	x			ALP
ALT		x	x	x					x	x	x			ALT
AMYL			x											AMYL
AST	x					x		x			x			AST
BUN		x	x	x		x	x	x	x	x	x			BUN
CA ⁺⁺	x	x	x			x	x	x			x			CA ⁺⁺
CHOL ⁻									x			x		CHOL ⁻
CHW		x												CHW
CK	x					x		x						CK
Cl ⁻				x	x			x			x			Cl ⁻
CREA		x	x	x		x	x			x	x			CREA
Gallsyra	x								x					Gallsyra
GGT						x		x					x	GGT
GLOB	x	x	x			x		x			x			GLOB
GLU	x	x	x	x		x	x			x	x			GLU
K ⁺	x		x	x		x	x				x			K ⁺
MG ⁺⁺								x						MG ⁺⁺
NA ⁺	x		x	x		x	x				x			NA ⁺
PHB													x	PHB
PHOS	x	x	x					x	x					PHOS
T4												x		T4
TBIL		x	x			x			x		x		x	TBIL
tCO ₂				x	x	x	x				x			tCO ₂
TP	x	x	x			x		x		x	x			TP
UA	x													UA



VetScan VS2 NG

TRIOLAB

Tel: 031-81 72 00
www.triolab.se
info@triolab.se

Alla rotorer kan användas till VetScan VS2 och VS2 NG. Rotorer märkta med * kan användas till VetScan Classic