

▶ VETSCAN i-STAT CARTRIDGE TEST MENU

A Wide Range of Cartridges Performing Advanced Diagnostics


		CG4+ #600-9000	CG8+ #600-9001	G #600-9002	Glu Strips #600-9011	Crea #600-9003	E3+ #600-9004	6+ #600-9005	CHEM8+ #600-9008	EC8+ #600-9007	ACT Celite® #600-9006	cTnl #600-9009
Hematology	Hematocrit (Hct)		•				•	•	•	•		
	Hemoglobin (Hgb)		•				•	•	•	•		
Chemistry	Blood Urea Nitrogen (BUN)							•	•	•		
	Creatinine					•			•			
	Ionized Calcium (iCa)		•						•			
	Glucose (Glu)		•	•	•			•	•	•		
Electrolytes	Chloride (Cl)*						•	•	•	•		
	Sodium (Na)		•				•	•	•	•		
	Potassium (K)		•				•	•	•	•		
Acid Base	pH	•	•							•		
	PCO ₂	•	•							•		
	HCO ₃	•	•							•		
	TCO ₂	•	•						•	•		
	Anion Gap								•	•		
	Base Excess	•	•							•		
Blood Gas	PO ₂	•	•									
	sO ₂	•	•									
Specialty	Lactate	•										
	ACT (Celite)										•	
	Cardiac Troponin I (cTnl)											•

* Available only on the VetScan i-STAT 1 Analyzer


▶ THREE STEPS TO ACCURATE RESULTS



1) Fill the cartridge to the fill line.



2) Insert cartridge.



3) Read results in 2 minutes.

FEATURES

- Small sample size – 95 µl or less
- Most results in 2 minutes
- Accurate and reliable results
- Low maintenance
- Automatic calibration
- Automated quality control checks
- Ergonomic design, handheld versatility, battery operated
- Stores up to 5,000 tests
- Rapid infrared printing through a convenient mini-printer (optional)
- Integrated barcode reader
- On board glucose strip reader
- Features a Cardiac Troponin I Cartridge

BENEFITS

- Perfect for either the anemic kitten or fractious cat
- Neither your client or you have to wait for results
- Allows you to make quick and confident treatment decisions
- Enhances usability in every clinical situation
- Prior patient results available in the exam room or in the field
- Simple, easy, convenient
- Quick general or emergency assessment with one cartridge

